



PPL PATIENT SPECIFIC - REUSABLE /DISPOSABLE SLINGS Weight Ranges

CAUTION - NEVER EXCEED WEIGHT CAPACITY OF PATIENT LIFT!

WEIGHT RANGE FOR ARJO LIFTS	PATIENT HEIGHT RANGE	BINDING COLOR	
PPLARJ10-S SMALL 80-130 LBS	4'10" TO 5'3"	RED	
PPLARJ10-M MEDIUM 120-180 LBS	5'3" TO 5'9"	YELLOW	
PPLARJ10-L LARGE 160-270 LBS	5'8" TO 6'	GREEN	
PPLARJ10-X EXTRA LARGE 250-360 LBS	5'11" TO 6'2"	PURPLE	
PPLARJ10-XX EXTRA EXTRA LARGE 350-420 LBS	6' TO 6'4"	WHITE	

WEIGHT RANGE FOR LOOP LIFTERS	PATIENT HEIGHT RANGE	BINDING COLOR	
PPL15-S SMALL 90-140 LBS	4'10" TO 5'3"	RED	
PPL15-M MEDIUM 130-180 LBS	5'4" TO 5'10"	YELLOW	
PPL15-L LARGE 170-280 LBS	5'8" TO 6'	GREEN	
PPL15-X EXTRA LARGE 270-410 LBS	5'10" TO 6'4"	PURPLE	
PPL15-XX EXTRAEXTRA LARGE 400-700LBS PLUS *SWL - SAFE WORKING LOAD 1000 LBS	6'2" TO 6'7"	WHITE	

SLING SIZING FOR PPL SLINGS WITH KEYHOLE CLIPS AND LOOPS

1. MEASURE PATIENT IN A SEATED POSITION. THE MEASUREMENT SHOULD BE FROM THE COCCYX TO THE TOP OF THE HEAD.

2. BODY PROPORTIONS AND PATIENT WIDTH SHOULD BE CONSIDERED BY CLINICAL STAFF BEFORE ASSIGNING SLING SIZE!