

### Beka Motion Arjo Parker Bath 500 Comparison

#### BEKA

- 1). Length (Entry/Tilting) 1.76mm / 1.908mm to **64 1/4" 75 1/8" Less Space**
- 2). Tub Width 860mm to **33 7/8" 4" wider**
- 3). Lowest entry ht 490mm to **19 1/4" 2" Lower**
- 4). Working Ht 810mm/1.010 to **31 3/4" to 40" 8" HT adj**
- 5). Water Drainage Liter/Min 75L to **20Gal**
- 6). Patient Weight 210kg to **463lbs**
- 7). Filling Volume (w/o Resident/W Resident) 270L/185L to **71 Gal/48 Gal**
- 8). Prefilling volume (Foot Zone) 98L to **26 Gal**
- 9). 10 Air Spa Jets

#### Parker

- 73" /77"**
- 29 1/2
- 21 1/2**
- 5" ht adj**
- 463lbs\***
- 26 Gal**
- 9 Spa Jets

#### Observations:

- Deeper Bath more Immersion
- Dual Handles internally for Additional security
- Dual control Handset and Panel (All controls)
- Can accommodate existing Maxi Move with DPS
- Door Lock System
- Lower seat height

- Single
- Hand set UP/Down Panel Water control
- NO
- NO Manual Latch

\*Parker Bath 500- 25.4 Gallons foot well water weight 203.2  
 Parker Bath 500 Maximum weight water and Patient 661lbs  
 Parker bath Max Patient Weight 457lbs

**UNITY-HEALTHCARE SOLUTIONS**  
**25967 CONIFER ROAD, CONIFER, CO 80421**  
**PH: 303-838-7320 ~ Fax: 303-838-4550**  
<http://www.unity-healthcare.com/>

