

UNITY-HEALTHCARE SOLUTIONS

Unity Healthcare provides your Disinfection Solutions

Complete Area Disinfection

With the air-assisted electrostatic spray-gun, applying chemicals efficiently and effectively is easy. The electrostatic charge ensures that even the hard to reach places are coated evenly. The sanitation line is the ideal solution for bio-security, fighting MRSA, E. Coli, salmonella, swine flu, C. Diff, H1N1, and various other bacteria and viruses.

SC-ET-HD-Electrostatic Sprayer



- Heavy-Duty, with an air-to-air heat exchanger for extended use
- Durable wheeled case with extendable handle
- Includes 1 MaxCharge™ Electrostatic Spray-gun
- Non-removable tank with drain hose and valve
- Self-contained, does not require an external air supply
- Standard 15 ft hose, also the maximum hose length for this system
- Internal tank holds enough mix for 1 hour
- Runs on 110 volts, 220-volt 50/60 Hertz available
- Compatible with all conventional, water-based chemicals and fungicides
- Quiet operation
- Spray range up to 8

SC-EB-Electrostatic Sprayer



- Compact, unobtrusive sprayer ideal for light-duty quick applications
- Polyresin case has wheels and extendable handle for easy transport
- Includes 1 MaxCharge™ electrostatic spray gun
- Oversized, semi-pneumatic tires
- 2 one-liter bottles for easy spot spraying
- Standard 30 ft hose, up to 75 ft maximum available
- Compatible with most conventional chemicals
- Quiet operation
- Runs on 110 volts, 220-volt 50/60 Hertz available
- Spray range up to 8 ft

SNB-Electrostatic Sprayer



- If you have your own compressor, SNB is a quick and less expensive option to spray electrostatically.
- This sprayer is ideal for hospitals, schools, hazmat locations, and fire department applications.
- Two one-liter Nalgene bottles are included.
- Bottles can be attached on a hip belt, or hand-held, for 15 minutes of spraying time each.
- Spraywand and Bottle (**compressor not included**)
- 45 feet of additional hose can be added to the standard hose length, providing a maximum length of 75 feet.
- Stainless steel connections (optional).

AUVS- Disinfectant Device



UV- Box Disinfectant Device lets you put the solution at the source with KR615 Germicidal Enclosure



- Decontaminate mobile hand held device surfaces, such as phones, tablets, and other non-critical devices.
- Working continuously between terminal cleanings and providing a better option than wipes.
- Patented UV technology delivers a kill rate greater than or equal to 99.99% for the toughest contaminants including MRSA and C. diff. spores.
- It protect staff, visitors and patients from cross contamination.
- Fills the gaps in your Infection Prevention Protocol.
- Working continuously between terminal cleanings and providing a better option than wipes•
- Eliminates cross contamination from handheld devices to hands
- Improves hand hygiene programs
- Easy to Use: One button operation
- Inexpensive
- Does not damage electronics
- Consistent results, cheaper and faster than disposable wipes
- High usage by staff eliminates contact dermatitis from disinfectants
- No waste or disposables

DISINFECT:



- Cell Phones
- Tablets
- Employee Badges
- Stethoscope
- TV Remote Controls and Hand Held Devices
- Project Folders

Unity-Healthcare Solutions, LLC
25967 Conifer Road; Conifer, CO 80433
Phone: 303-838-7320
Fax: 303-838-4550



NoroxyCDIFF:

- EPA Registered Broad Spectrum

Efficacy:

C.diff, MRSA, VREHBV, Salmonella, HIV, Influenza A: H1N1 E.Coli, Norovirus, Klebsiella Pneumonia, Pseudomonas, Aeruginosa, Acinetobacter Baumannii, Multi-drug resistant bacteria, Other bacteria, viruses and fungi.

- Disinfectant without unwanted side effects
- FDA Clearance: As Sanitizers on Food Contact Surfaces Bactericide and Viricidal
- One step Hospital use sporicidal Disinfectant Cleaner
- Formatted to prevent cross contamination on Hard Non-Porous Surfaces
- Improved room turn over times

Can be applied using STS cleaning

Made in U.S.A. •EPA Registered

Broad Spectrum Microbial Action

FDA Clearance as Sanitizer on Food Contact Surfaces

Cdiff Kill Time 2 Minutes (Most Difficult)

Active Ingredients:

Hydrogen Peroxide, Peroxyacetic Acid

Outstanding Disinfecting Properties

Effective Against:

Bacteria Viruses, Spores and Fungi

1 of only 50 EPA Registered Products Effective Against CDiff42 are Formulated with Same Form of Bleach/Chlorine

Remaining Products Registered are Either Peracetic Acid Liquid or Hydrogen Peroxide •NOROxyCdiff Penetrates cell Membrane to Disrupt Enzyme to Kill Microorganism Protocols

Effective against Bacteria, Viruses, Spores and Fungus

Unique ability to penetrate Acetic acid Irreversibly

Disrupts Enzyme System and the Destruction of the microorganism.

NO Bleach, NO Ammonia, NO Phosphates- Non-Corrosive- No rinse, No wipe



MediDefense Penetrex

- Biostatic surface treatment reduces risks of germ-transfer
- Around the clock germ protection
- Environmentally safe, non-toxic, and non-leaching
- Supports patient protection mandates



Disinfectant – 275

- EPA registered
- Non-toxic,
- Biodegradable,
- Safe, fast-acting healthcare-grade cleaners, and disinfectants.
- These electrochemical activation solutions clean and degrease surfaces
- Destroy pathogens, viruses, algae, fungi
- Non-toxic, eco-friendly, fast acting and safe to use
- Ready to use Hospital Disinfectant
- Fast acting, residue-free, non-toxic disinfectant
- Effective against Pseudomonas aeruginosa, Salmonella enterica, Staphylococcus aureus, Swine influenza Virus (H1N1)



Disinfectant – 525

- EPA registered
- Non-toxic,
- Biodegradable,
- Safe, fast-acting healthcare-grade cleaners, and disinfectants.
- These electrochemical activation solutions clean and degrease surfaces
- Destroy pathogens, viruses, algae, fungi
- Non-toxic, eco-friendly, fast acting and safe to use
- Ready to use Broad Spectrum Hospital Disinfectant
- Effective against Methicillin Resistant, Staphylococcus Aureus (MRSA), Swine Influenza Virus (H1V1), Escherichia coli (NDM), Vancomycin resistant enterococci (VRE), Human Immunodeficiency Virus Type 1 (HIV-1), Mycobacterium bovis (BCG), Salmonella enterica
- Hypochlorous acid (HOCL) kills CRE's on contact
- Fast-acting, residue free, non-toxic disinfectant.
- Micro-organism cannot develop resistance against Hypochlorous acid (HOCL), unlike the case with chemicals (e.g. occurrence of super bugs)

mPulse HAND SANITIZER

Eliminates 99.99% of germs on contact

Active Antimicrobial Protection increases germ killing effectiveness with continued use.

- Active ingredient is categorized as a quaternary ammonium and meets the EPA's criteria for use against SARS-CoV-2, the virus that causes COVID-19
- Extends protection between hand washings
- Alcohol free moisturizing formula will not crack your hands or strip away your skin's natural oils
- Protects against dirt, oil and grease build up between hand washings
- Approved for food handling without the need to rinse



mPulse Hand Sanitizer is the perfect companion to your KR615 "UV Box"

Unity-Healthcare Solutions, LLC
25967 Conifer Road; Conifer, CO 80433
Phone: 303-838-7320
Fax: 303-838-4550