



**PATIENT SPECIFIC - REUSABLE / DISPOSABLE
 SLINGS Weight Ranges**

CAUTION - NEVER EXCEED WEIGHT CAPACITY OF PATIENT LIFT!

| WEIGHT RANGE FOR CLIP STYLE SLINGS | PATIENT HEIGHT RANGE | BINDING COLOR | |
|---|-----------------------------|----------------------|---|
| SMALL 80-130 LBS | 4'10" TO 5'3" | RED |  |
| MEDIUM 120-180 LBS | 5'3" TO 5'9" | YELLOW | |
| LARGE 160-270 LBS | 5'8" TO 6'0" | GREEN | |
| EXTRA LARGE 250-360 LBS | 5'11" TO 6'2" | PURPLE | |
| EXTRA EXTRA LARGE 350- 500 LBS | 6'0" TO 6'4" | WHITE | |

| WEIGHT RANGE FOR LOOP SLINGS | PATIENT HEIGHT RANGE | BINDING COLOR | |
|-------------------------------------|-----------------------------|----------------------|---|
| SMALL 90-140 LB | 4'10" TO 5'3" | RED |  |
| MEDIUM 130-180 LBS | 5'4" TO 5'10" | YELLOW | |
| LARGE 170-280 LBS | 5'8" TO 6'0" | GREEN | |
| EXTRA LARGE 270-410 LBS | 5'10" TO 6'4" | PURPLE | |
| *EXTRA EXTRA LARGE 400-700LBS PLUS | 6'2" TO 6'7" | WHITE | |

*SWL - SAFE WORKING LOAD 1000 LBS

SLING SIZING FOR PPL SLINGS WITH KEYHOLE CLIPS AND LOOPS

1. MEASURE PATIENT IN A SEATED POSITION. THE MEASUREMENT SHOULD BE FROM THE COCCYX TO THE TOP OF THE HEAD.
2. BODY PROPORTIONS AND PATIENT WIDTH SHOULD BE CONSIDERED BY CLINICAL STAFF BEFORE ASSIGNING SLING SIZE!

Unity-Healthcare, LCC
25967 Conifer Road; Conifer, CO 80433
Main Line: 303-838-7320 Fax: 303-838-4550

